



Product description

| | |
|---------------------|--|
| Order number | 83.3944 |
| Product description | Pressure filtration unit, Filtropur L 0.2 S, pre-filter: GF, membrane: PES, filtration surface: 20 cm ² , pore size: 0.2 µm, for sterile filtration, connection: Tube connector, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag |

Product characteristics

| | |
|------------------------|------------------------|
| Example of use | for sterile filtration |
| Filter | yes |
| Surface area of filter | 20 cm ² |
| Pore size | 0,2 |
| Connection | Tube connector |

Size

| | |
|-------------------|----------|
| Filtration volume | 10000 ml |
|-------------------|----------|

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|------------------------|
| Product material | Polycarbonate (PC) |
| Membrane material | Polyethersulfone (PES) |

Purity & certification

| | |
|------------------|--|
| Product category | no medical device IVD |
| Purity standard | sterile, non-pyrogenic/endotoxin-free, non-cytotoxic |
| Sterilisation | Gas sterilisation |
| Non-cytotoxic | yes |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 10 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 10 |
| Piece (s) / outer case | 50 |
| Piece(s) / pallet | 400 |
| Depth of box | 352 mm |
| Width of box | 155 mm |
| Height of box | 123 mm |
| Depth of case | 530 mm |
| Width of case | 370 mm |
| Height of case | 280 mm |
| Case volume | 0.0549 cbm |
| Weight of product | 0.025 kg |
| Weight of case | 4.24 kg |
| EAN of inner box | 4038917926225 |
| EAN case | 4038917926232 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

