



Product description

| | |
|---------------------|---|
| Order number | 94.6077.305 |
| Product description | lumox® dish 50, Tissue culture dish, with foil base, Ø: 50 mm, for suspension cells, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 5 piece(s)/bag |

Product characteristics

| | |
|----------------|-----------------------|
| Closure type | slip-on lid |
| Design | with foil base |
| Cell type | suspension cells |
| Base shape | flat base |
| Growth area | 20.35 cm ² |
| Volume of work | 5 ml - 10 ml |
| Vessel type | Dish |
| Coding colour | green |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
lumox® dish 50, Tissue culture dish, with foil base, Ø: 50 mm, suspension cells

Page 2

Size

| | |
|-------------------|-------|
| Diameter | 50 mm |
| Height of product | 12 mm |

Material & colours

| | |
|------------------|------------------|
| Product material | lumox® |
| Frame material | Polystyrene (PS) |

Purity & certification

| | |
|------------------|--|
| Product category | no medical device IVD |
| Purity standard | sterile, non-pyrogenic/endotoxin-free, non-cytotoxic |
| Sterilisation | Electron irradiation |
| Pyrogen-free | yes |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 50 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 200 |
| Piece(s) / pallet | 7200 |
| Depth of case | 362 mm |
| Width of case | 162 mm |
| Height of case | 125 mm |
| Case volume | 0.0073 cbm |
| Weight of product | 0.01 kg |
| Weight of case | 3.52 kg |
| EAN of inner box | 4038917256070 |
| EAN case | 4038917256063 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT