




Product Specification
lumox® dish 35, Tissue culture dish, with foil base, Ø: 35 mm, Adherent cells

Page 1



Product description

Order number	94.6077.331
Product description	lumox® dish 35, Tissue culture dish, with foil base, Ø: 35 mm, for Adherent cells, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 5 piece(s)/bag

Product characteristics

Closure type	slip-on lid
Design	with foil base
Cell type	Adherent cells
Base shape	flat base
Growth area	6.35 cm ²
Volume of work	2,5 ml
Vessel type	Dish
Coding colour	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
lumox® dish 35, Tissue culture dish, with foil base, Ø: 35 mm, Adherent cells

Page 2

Size

Diameter	35 mm
Height of product	6 mm

Material & colours

Product material	lumox®
Frame material	Polystyrene (PS)

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	50
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	250
Piece(s) / pallet	15000
Depth of case	348 mm
Width of case	253 mm
Height of case	265 mm
Case volume	0.0233 cbm
Weight of product	0.0037 kg
Weight of case	2.92 kg
EAN of inner box	4038917255110
EAN case	4038917255103

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT