



### **Product description**

Order number 94.6001.061  
Product description Universal turning device, for 4 miniPERM® bioreactors, not autoclavable, 115/230 V, 1 piece(s)/case

### **Product characteristics**

Suitable for 4 miniPERM® bioreactors

### **Size**

Width of product 26 cm  
Length of product 46 cm

### **Purity & certification**

Product category no medical device | IVD  
Batched no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	30
Depth of case	398 mm
Width of case	358 mm
Height of case	325 mm
Case volume	0.0463 cbm
Weight of product	5.861 kg
Weight of case	6.3 kg
EAN of inner box	4038917254878
EAN case	4038917254861

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

