



### **Product description**

Order number	94.6001.153
Product description	miniPERM® , Nutrient module, for miniPERM® production module, autoclavable up to 121 °C, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

### **Product characteristics**

Closure type	screw cap
Suitable for	miniPERM® production module
Volume of work	350 ml - 450 ml
Autoclavable (max. °C)	121 °C

### **Size**

Sample volume	400 ml
Diameter	102 mm
Height of product	129 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polycarbonate (PC)
Closure material	High Density Polyethylene (HD-PE)

### **Purity & certification**

Product category	no medical device   IVD
Purity standard	non-pyrogenic/endotoxin-free, non-cytotoxic
Batched	yes

### **Packaging**

Minimum order qty.	4
Type of smallest subpackaging	bag
Piece(s) / inner box	4
Piece (s) / outer case	24
Piece(s) / pallet	480
Depth of case	498 mm
Width of case	364 mm
Height of case	310 mm
Case volume	0.0562 cbm
Weight of product	0.1208 kg
Weight of case	4.74 kg
EAN of inner box	4038917254991
EAN case	4038917254984

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

