




Product Specification
miniPERM® , Disposable syringe 2 ml, for
miniPERM® bioreactor, Luer

Page 1



Product description

Order number	94.6077.136
Product description	miniPERM® , Disposable syringe 2 ml, for miniPERM® bioreactor, Luer, not autoclavable, sterile, 1 piece(s)/blister

Product characteristics

Suitable for	miniPERM® bioreactor
Connection	Luer

Size

Max. volume	2 ml
-------------	------

Material & colours

Product material	Polypropylene (PP)
------------------	--------------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
miniPERM® , Disposable syringe 2 ml, for
miniPERM® bioreactor, Luer

Page 2

Purity & certification

Product category	Medical device
CE certified	CE supplier
Purity standard	sterile
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	100
Type of smallest subpackaging	blister
Piece(s) / inner box	100
Piece (s) / outer case	100
Piece(s) / pallet	32400
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.003 kg
Weight of case	0.42 kg
EAN of inner box	4038917255929
EAN case	4038917255912

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT