



### **Product description**

|                     |   |
|---------------------|---|
| Order number        | 94.6140.402   |
| Product description | x-well cell culture chamber, 4 wells, on PCA slide, removable frame, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 6 piece(s)/blister |

### **Product characteristics**

|                 |                     |
|-----------------|---------------------|
| Closure type    | slip-on lid         |
| Design          | on PCA slide        |
| Cell type       | Adherent cells      |
| Number of wells | 4                   |
| Growth area     | 1.9 cm <sup>2</sup> |
| Volume of work  | 1 ml                |
| Removable frame | yes                 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

|                   |       |
|-------------------|-------|
| Width of product  | 24 mm |
| Length of product | 76 mm |

### **Material & colours**

|                  |                        |
|------------------|------------------------|
| Product material | Polycycloalkanes (PCA) |
| Frame material   | Polystyrene (PS)       |

### **Purity & certification**

|                  |  |
|------------------|--|
| Product category | no medical device   IVD                              |
| Purity standard  | sterile, non-pyrogenic/endotoxin-free, non-cytotoxic |
| Sterilisation    | Gas sterilisation                                    |
| Pyrogen-free     | yes  |
| Batched          | yes  |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 6             |
| Type of smallest subpackaging | blister       |
| Piece(s) / inner box          | 6             |
| Piece (s) / outer case        | 96            |
| Piece(s) / pallet             | 10368         |
| Depth of box                  | 127 mm        |
| Width of box                  | 19 mm         |
| Height of box                 | 258 mm        |
| Depth of case                 | 275 mm        |
| Width of case                 | 270 mm        |
| Height of case                | 185 mm        |
| Case volume                   | 0.0137 cbm    |
| Weight of product             | 0.0143 kg     |
| Weight of case                | 2.25 kg       |
| EAN of inner box              | 4038917303293 |
| EAN case                      | 4038917303286 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

