



Product description

Order number	94.6190.102
Product description	x-well cell culture chamber, 1 well, on coverglass, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 6 piece(s)/blister

Product characteristics

Closure type	slip-on lid
Design	on coverglass
Cell type	Adherent cells
Number of wells	1
Growth area	9 cm ²
Volume of work	4 ml
Removable frame	no

Size

Width of product	24 mm
Length of product	60 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material borosilicate glass
Frame material Polystyrene (PS)

Purity & certification

Product category no medical device | IVD
Purity standard sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation Gas sterilisation
Pyrogen-free yes
Batched yes

Packaging

Minimum order qty. 6
Type of smallest
subpackaging blister
Piece(s) / inner box 6
Piece (s) / outer case 96
Piece(s) / pallet 10368
Depth of box 127 mm
Width of box 19 mm
Height of box 258 mm
Depth of case 275 mm
Width of case 270 mm
Height of case 185 mm
Case volume 0.0137 cbm
Weight of product 0.009 kg
Weight of case 1.74 kg
EAN of inner box 4038917256674
EAN case 4038917256667

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

