



### Product description

Order number	82.1972.700
Product description	Deep Well plate, Low DNA-binding, 96 well, 2.2 ml, conical wedge, without label, PCR Performance Tested, material: PP, Edged cavities, 4 piece(s)/bag

### Product characteristics

Label/ Print	without label
Closure type	film
Application	collection and storage of diagnostic specimens
Base shape	conical wedge
Number of wells	96
Volume of work	1 $\mu$ l - 2,2 ml
Centrifugation max (RCF)	4500 x g
Centrifugation temperature	22
Centrifugation time	5 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### Size

Sample volume	2.2 ml
Max. volume	2.2 ml
Diameter	7.52 mm
Width of product	85 mm
Height of product	41.85 mm
Length of product	127.76 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

### Purity & certification

Product category	no medical device   IVD
Purity standard	PCR Performance Tested
Batched	yes

### Packaging

Minimum order qty.	24
Type of smallest subpackaging	bag
Piece(s) / inner box	4
Piece (s) / outer case	24
Piece(s) / pallet	2304
Depth of case	283 mm
Width of case	268 mm
Height of case	210 mm
Case volume	0.0159 cbm
Weight of product	0.0868 kg
Weight of case	2.43 kg
EAN of inner box	4038917958042
EAN case	4038917958059

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

