



### **Product description**

Order number	95.2093
Product description	Film, free from DNase/RNase, transparent, pierceable, application area: -80 - 120 °C

### **Product characteristics**

Closure type	film
Light protection	no
Suitable for qPCR	no
Minimum temperature	-80 °C
Maximum temperature	120 °C

### **Size**

Width of product	75 mm
Length of product	118 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Film, free from DNase/RNase, transparent,  
pierceable**

**Page 2**

**Material & colours**

Colour of product                      transparent

**Purity & certification**

Product category                      no medical device | IVD  
Purity standard                        free from DNase/RNase  
Batched                                    yes

**Packaging**

Minimum order qty.                    20  
Piece(s) / inner box                    20  
Piece (s) / outer case                  240  
Piece(s) / pallet                        240  
Case volume                              0.0142 cbm  
Weight of product                      0.012 kg  
Weight of case                          3 kg  
EAN of inner box                        4038917335706  
EAN case                                  4038917335690

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**