



Product description

Order number	86.1203
Product description	Disposal bags, volume: 40 l, (LxW): 780 x 600 mm, material: PP, transparent, film thickness: 50 µm, autoclavable up to 134 °C, with print, print Biohazard, label/print: red, 50 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	red
Liquid-tight	yes
Autoclavable (max. °C)	134 °C

Size

Width of product	600 mm
Length of product	780 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Disposal bags, 40 l, (LxW): 780 x 600 mm, PP,
transparent, with print Biohazard

Page 2

Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Latex-free yes
Safety data sheet no
Batched no

Packaging

Minimum order qty. 250
Type of smallest bag
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 250
Piece(s) / pallet 6000
Depth of case 290 mm
Width of case 110 mm
Height of case 630 mm
Case volume 0.0201 cbm
Weight of product 0.0415 kg
Weight of case 10.76 kg
EAN of inner box 4038917085779
EAN case 4038917004404

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT