



## Product description

Order number	86.1203.103
Product description	Disposal bags, volume: 40 l, (LxW): 780 x 600 mm, material: PP, yellow, film thickness: 50 µm, autoclavable up to 134 °C, with print, label/print: blue, 50 piece(s)/bag

## Product characteristics

Label/ Print	with print
Colour of print/label	blue
Liquid-tight	yes
Autoclavable (max. °C)	134 °C

## Size

Width of product	600 mm
Length of product	780 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Disposal bags, 40 l, (LxW): 780 x 600 mm, PP,  
yellow, with print**

**Page 2**

**Material & colours**

Product material Polypropylene (PP)  
Colour of product yellow

**Purity & certification**

Product category no medical device | IVD  
Latex-free yes  
Safety data sheet no  
Batched no

**Packaging**

Minimum order qty. 250  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 50  
Piece (s) / outer case 250  
Piece(s) / pallet 7000  
Depth of case 623 mm  
Width of case 288 mm  
Height of case 115 mm  
Case volume 0.0206 cbm  
Weight of product 0.0415 kg  
Weight of case 10.7 kg  
EAN of inner box 4038917126649  
EAN case 4038917020732

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**