



## Product description

Order number	86.1206.103
Product description	Disposal bags, volume: 80 l, (LxW): 1,120 x 700 mm, material: PP, yellow, film thickness: 50 µm, autoclavable up to 134 °C, with print, print Biohazard, label/print: blue, 50 piece(s)/bag

## Product characteristics

Label/ Print	with print
Colour of print/label	blue
Liquid-tight	yes
Autoclavable (max. °C)	134 °C

## Size

Width of product	700 mm
Length of product	1.120.00 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Disposal bags, 80 l, (LxW): 1,120 x 700 mm, PP,  
yellow, with print Biohazard**

**Page 2**

**Material & colours**

Product material Polypropylene (PP)  
Colour of product yellow

**Purity & certification**

Product category no medical device | IVD  
Latex-free yes  
Safety data sheet no  
Batched no

**Packaging**

Minimum order qty. 150  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 50  
Piece (s) / outer case 150  
Piece(s) / pallet 3600  
Depth of case 378 mm  
Width of case 122 mm  
Height of case 810 mm  
Case volume 0.0374 cbm  
Weight of product 0.072 kg  
Weight of case 11.38 kg  
EAN of inner box 4038917160452  
EAN case 4038917160445

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-06-30

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**