



Product description

Order number	82.1581.001
Product description	Micro test plate, 96 well, slip-on lid, flat base, lid assembled, material: PS, transparent, free from DNA/DNase/RNase, free from pyrogens/endotoxins, non-cytotoxic, sterile, 1 piece(s)/blister

Product characteristics

Label/ Print	with plastic label
Closure type	slip-on lid
Base shape	flat base
Number of wells	96
Vessel type	Plate

Size

Max. volume	390 μ l
Width of product	85.2 mm
Height of product	16.55 mm
Length of product	127.8 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard free from DNA/DNase/RNase, free from
pyrogens/endotoxins, non-cytotoxic, sterile
Sterilisation Gamma irradiation
Batched yes

Packaging

Minimum order qty. 50
Type of smallest blister
subpackaging
Piece(s) / inner box 1
Piece (s) / outer case 50
Piece(s) / pallet 3200
Depth of case 332 mm
Width of case 297 mm
Height of case 215 mm
Case volume 0.0212 cbm
Weight of product 0.0568 kg
Weight of case 3.54 kg
EAN of inner box 4038917131056
EAN case 4038917347952

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

