



Product description

Order number	82.1581.220
Product description	ELISA plate, 96 well, flat base, without lid, material: PS, black, High Binding, free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic, 25 piece(s)/bag

Product characteristics

Label/ Print	with print
Base shape	flat base
Number of wells	96
Vessel type	Plate

Size

Max. volume	390 µl
Width of product	85.2 mm
Height of product	14.35 mm
Length of product	127.8 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
ELISA plate, 96 well, flat base, PS, black, High Binding

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product black

Purity & certification

Product category no medical device | IVD
Purity standard free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic
Pyrogen-free yes
Batched yes

Packaging

Minimum order qty. 50
Type of smallest subpackaging bag
Piece(s) / inner box 25
Piece (s) / outer case 50
Piece(s) / pallet 2000
Depth of case 279 mm
Width of case 232 mm
Height of case 365 mm
Case volume 0.0236 cbm
Weight of product 0.0408 kg
Weight of case 2.86 kg
EAN of inner box 4038917365390
EAN case 4038917365383

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT