



### **Product description**

Order number 71.9923.212  
Product description Weigh tray, 70 ml, (LxW): 72 x 72 mm, material: PVC, white, 100 piece(s)/bag

### **Size**

Sample volume 70 ml  
Max. volume 100 ml  
Width of product 72 mm  
Length of product 72 mm

### **Material & colours**

Product material Polyvinylchloride (PVC)  
Colour of product white

### **Purity & certification**

Product category no medical device | IVD  
Batched yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	2000
Piece(s) / pallet	70000
Depth of case	428 mm
Width of case	234 mm
Height of case	358 mm
Case volume	0.0359 cbm
Weight of product	0.0024 kg
Weight of case	5.6 kg
EAN of inner box	4038917153669
EAN case	4038917038249

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

