



### **Product description**

Order number	84.1027.010
Product description	Adhesive label, (LxW): 62 x 40 mm, material: paper, white, label layout: Name + date, print: white/black, 5,000 piece(s)/roll

### **Product characteristics**

Colour of print/label	white/black
-----------------------	-------------

### **Size**

Width of product	40 mm
Length of product	62 mm

### **Material & colours**

Product material	paper
Colour of product	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

### **Packaging**

Minimum order qty.                      5000  
Type of smallest  
subpackaging                              roll  
Piece(s) / inner box                      5000  
Piece (s) / outer case                      20000  
Piece(s) / pallet                            1280000  
Depth of case                              398 mm  
Width of case                                298 mm  
Height of case                               215 mm  
Case volume                                0.0255 cbm  
Weight of product                        0.0002 kg  
Weight of case                              4.44 kg  
EAN of inner box                         4038917124706  
EAN case                                      4038917019354

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

