



Product description

Order number	84.1556.010
Product description	Adhesive label, (LxW): 40 x 26 mm, material: paper, white, label layout: Name + date, print: white/black, 5,000 piece(s)/roll

Product characteristics

Colour of print/label	white/black
-----------------------	-------------

Size

Width of product	26 mm
Length of product	40 mm

Material & colours

Product material	paper
Colour of product	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 5000
Type of smallest
subpackaging roll
Piece(s) / inner box 5000
Piece (s) / outer case 30000
Piece(s) / pallet 2880000
Depth of case 333 mm
Width of case 303 mm
Height of case 150 mm
Case volume 0.0151 cbm
Weight of product 0.0001 kg
Weight of case 2.85 kg
EAN of inner box 4038917130578
EAN case 4038917023139

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

