



Product description

Order number 84.2620.003
Product description Adhesive label, (LxW): 30 x 20 mm, 3,500 piece(s)/roll

Size

Width of product 20 mm
Length of product 30 mm

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 7000
Type of smallest subpackaging roll
Piece(s) / inner box 3500
Piece (s) / outer case 7000
Piece(s) / pallet 2268000
Depth of case 503 mm
Width of case 375 mm
Height of case 155 mm
Case volume 0.0292 cbm
Weight of product 0.0001 kg
Weight of case 0.73 kg
EAN of inner box 4038917900362
EAN case 4038917900379

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

