



Product description

Order number 95.1760.100
Product description Rubber foot, for pipette tip boxes, grey

Product characteristics

Suitable for pipette tip boxes

Material & colours

Product material Thermoplastic elastomers (TPE)
Colour of product grey

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 25
Piece(s) / inner box 25
Piece (s) / outer case 2000
Piece(s) / pallet 648000
Depth of case 268 mm
Width of case 191 mm
Height of case 100 mm
Case volume 0.0051 cbm
Weight of product 0.0001 kg
Weight of case 0.45 kg
EAN of inner box 4038917522670
EAN case 4038917522663

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

