



Product description

Order number	70.1183.002
Product description	Pipette tip, working volume: 5 ml, transparent, filling level rings, suitable for Eppendorf, 250 piece(s)/bag

Product characteristics

Graduation	yes
Suitable for	Eppendorf
Volume of work	5 ml
Filter	no
Autoclavable (max. °C)	121 °C

Size

Max. volume	5 ml
Length of product	123 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 2500
Type of smallest subpackaging bag
Piece(s) / inner box 250
Piece (s) / outer case 2500
Piece(s) / pallet 25000
Depth of case 298 mm
Width of case 464 mm
Height of case 698 mm
Case volume 0.0966 cbm
Weight of product 0.0026 kg
Weight of case 7.3 kg
EAN of inner box 4038917126458
EAN case 4038917574594

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

