



Product description

Order number	70.1186
Product description	Pipette tip, working volume: 1,250 µl, transparent, filling level rings, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design, 384 piece(s)/bag

Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design
Volume of work	1.250 µl
Filter	no
Coding colour	blue
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

Max. volume 1250 µl
Length of product 103 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 3840
Type of smallest subpackaging bag
Piece(s) / inner box 384
Piece (s) / outer case 3840
Piece(s) / pallet 134400
Depth of case 347 mm
Width of case 258 mm
Height of case 443 mm
Case volume 0.0397 cbm
Weight of product 0.0011 kg
Weight of case 4.57 kg
EAN of inner box 4038917225236
EAN case 4038917225229

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

