



Product description

Order number	70.1187.001
Product description	Pipette tip, working volume: 10 ml, transparent, filling level rings, suitable for Eppendorf, Gilson and Brand and products of identical design, 100 piece(s)/bag

Product characteristics

Graduation	yes
Suitable for	Eppendorf, Gilson and Brand and products of identical design
Volume of work	10 ml
Filter	no
Autoclavable (max. °C)	121 °C

Size

Max. volume	10 ml
Length of product	160.8 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest subpackaging bag
Piece(s) / inner box 100
Piece (s) / outer case 1000
Piece(s) / pallet 20000
Depth of case 467 mm
Width of case 289 mm
Height of case 387 mm
Case volume 0.0522 cbm
Weight of product 0.0038 kg
Weight of case 6.6 kg
EAN of inner box 4038917574617
EAN case 4038917574600

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

