



Product description

Order number	70.1187.102
Product description	Pipette tip, working volume: 10 ml, blue, filling level rings, suitable for Sarpette® M, 25 piece(s)/box

Product characteristics

Graduation	yes
Suitable for	Sarpette® M
Volume of work	10 ml
Filter	no
Autoclavable (max. °C)	121 °C

Size

Max. volume	10 ml
Length of product	144 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product blue

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 250
Type of smallest subpackaging box
Piece(s) / inner box 25
Piece (s) / outer case 250
Piece(s) / pallet 7500
Depth of case 603 mm
Width of case 303 mm
Height of case 152 mm
Case volume 0.0278 cbm
Weight of product 0.0128 kg
Weight of case 3.6 kg
EAN of inner box 4038917951494
EAN case 4038917951500

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

