



Product description

Order number	70.1189.105
Product description	Pipette tip, working volume: 200 µl, transparent, filling level rings, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpiquette, Biohit and products of identical design, extra long, 96 piece(s)/box

Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpiquette, Biohit and products of identical design
Volume of work	200 µl
Filter	no
Coding colour	yellow
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Pipette tip, extra long, 200 µl, transparent, 96
piece(s)/box**

Page 2

Size

Max. volume 200 µl
Length of product 77 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 960
Type of smallest box
subpackaging
Piece(s) / inner box 96
Piece (s) / outer case 960
Piece(s) / pallet 120960
Depth of case 263 mm
Width of case 121 mm
Height of case 413 mm
Case volume 0.0131 cbm
Weight of product 0.0015 kg
Weight of case 1.94 kg
EAN of inner box 4038917314497
EAN case 4038917314480

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT