



### **Product description**

Order number	70.1189.215
Product description	Filter tip, working volume: 200 µl, transparent, filling level rings, Biosphere® , suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson and Finnpiquette and products of identical design, 96 piece(s)/box

### **Product characteristics**

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson and Finnpiquette and products of identical design
Volume of work	200 µl
Filter	yes
Coding colour	yellow

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Filter tip, 200 µl, transparent, Biosphere® , 96  
piece(s)/box**

**Page 2**

**Size**

Max. volume 200 µl  
Length of product 77 mm

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent  
Filter material Polyethylene (PE)

**Purity & certification**

Product category no medical device | IVD  
Purity standard Biosphere®  
Sterilisation Gamma irradiation  
Pyrogen-free yes  
Batched yes

**Packaging**

Minimum order qty. 960  
Type of smallest subpackaging box  
Piece(s) / inner box 96  
Piece (s) / outer case 960  
Piece(s) / pallet 120960  
Depth of case 241 mm  
Width of case 89 mm  
Height of case 406 mm  
Case volume 0.0087 cbm  
Weight of product 0.0014 kg  
Weight of case 1.94 kg  
EAN of inner box 4038917314589  
EAN case 4038917314572

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**