



Product description

Order number	70.1190.100
Product description	Gelloader pipette tip, working volume: 200 µl, transparent, suitable for common pipettors, 96 piece(s)/box

Product characteristics

Graduation	no
Example of use	Gelloader
Suitable for	common pipettors
Volume of work	200 µl
Filter	no
Coding colour	yellow
Autoclavable (max. °C)	121 °C

Size

Max. volume	200 µl
Length of product	67 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Gelloader pipette tip, 200 µl, transparent, 96
piece(s)/box**

Page 2

Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 960
Type of smallest box
subpackaging
Piece(s) / inner box 96
Piece (s) / outer case 960
Piece(s) / pallet 103680
Depth of case 475 mm
Width of case 255 mm
Height of case 100 mm
Case volume 0.0121 cbm
Weight of product 0.0015 kg
Weight of case 1.73 kg
EAN of inner box 4038917344593
EAN case 4038917344586

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT