



Product description

| | |
|---------------------|--|
| Order number | 70.3010.100 |
| Product description | Pipette tip, working volume: 10 µl, transparent, filling level rings, PCR Performance Tested, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design, 480 piece(s)/StackPack |

Product characteristics

| | |
|----------------|---|
| Graduation | yes |
| Suitable for | SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design |
| Volume of work | 10 µl |
| Filter | no |
| Coding colour | grey |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Pipette tip, 10 µl, transparent, PCR Performance
Tested, 480 piece(s)/StackPack

Page 2

Size

Max. volume 10 µl
Length of product 31 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard PCR Performance Tested
Batched yes

Packaging

Minimum order qty. 1920
Type of smallest StackPack
subpackaging
Piece(s) / inner box 1920
Piece (s) / outer case 1920
Piece(s) / pallet 307200
Depth of case 290 mm
Width of case 143 mm
Height of case 225 mm
Case volume 0.0093 cbm
Weight of product 0.0004 kg
Weight of case 1.69 kg
EAN of inner box 4038917894432
EAN case 4038917894425

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT