



### **Product description**

Order number	70.3030.020
Product description	Pipette tip, working volume: 200 µl, yellow, filling level rings, PCR Performance Tested, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 1,000 piece(s)/bag

### **Product characteristics**

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design
Volume of work	200 µl
Filter	no
Coding colour	yellow
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Pipette tip, 200 µl, yellow, PCR Performance**  
**Tested, 1,000 piece(s)/bag**

**Page 2**

**Size**

Max. volume 200 µl  
Length of product 51 mm

**Material & colours**

Product material Polypropylene (PP)  
Colour of product yellow

**Purity & certification**

Product category In vitrodiagnostic device, CE  
CE certified CE  
Purity standard PCR Performance Tested  
Batched yes

**Packaging**

Minimum order qty. 10000  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 1000  
Piece (s) / outer case 10000  
Piece(s) / pallet 480000  
Depth of case 388 mm  
Width of case 378 mm  
Height of case 215 mm  
Case volume 0.0315 cbm  
Weight of product 0.0003 kg  
Weight of case 3.83 kg  
EAN of inner box 4038917525985  
EAN case 4038917525978

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-01

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**