




**Product Specification**  
**Pipette tip, 200 µl, transparent, Biosphere® plus,  
96 piece(s)/box**

**Page 1**



**Product description**

Order number	70.3030.205
Product description	Pipette tip, working volume: 200 µl, transparent, filling level rings, Biosphere® plus, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 96 piece(s)/box

**Product characteristics**

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design
Volume of work	200 µl
Filter	no
Coding colour	yellow
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Pipette tip, 200 µl, transparent, Biosphere® plus,**  
**96 piece(s)/box**

**Page 2**

**Size**

Max. volume 200 µl  
Length of product 51 mm

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

**Purity & certification**

Product category In vitrodiagnostic device, CE  
CE certified CE  
Purity standard Biosphere® plus  
Sterilisation Gas sterilisation  
Batched yes

**Packaging**

Minimum order qty. 1920  
Type of smallest box  
subpackaging  
Piece(s) / inner box 480  
Piece (s) / outer case 1920  
Piece(s) / pallet 96000  
Depth of box 330 mm  
Width of box 111 mm  
Height of box 148 mm  
Depth of case 348 mm  
Width of case 325 mm  
Height of case 245 mm  
Case volume 0.0277 cbm  
Weight of product 0.0011 kg  
Weight of case 5.16 kg  
EAN of inner box 4038917847193  
EAN case 4038917847209

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**