



### **Product description**

Order number	70.3031.210
Product description	Pipette tip, working volume: 250 µl, transparent, filling level rings, PCR Performance Tested, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 96 piece(s)/box

### **Product characteristics**

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design
Volume of work	250 µl
Filter	no
Coding colour	yellow
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Pipette tip, 250 µl, transparent, PCR Performance Tested, Low retention, 96 piece(s)/box**

**Page 2**

**Size**

Max. volume	250 µl
Length of product	51 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

**Packaging**

Minimum order qty.	1920
Type of smallest subpackaging	box
Piece(s) / inner box	480
Piece (s) / outer case	1920
Piece(s) / pallet	96000
Depth of box	348 mm
Width of box	129 mm
Height of box	166 mm
Depth of case	348 mm
Width of case	325 mm
Height of case	245 mm
Case volume	0.0277 cbm
Weight of product	0.0011 kg
Weight of case	5.16 kg
EAN of inner box	4038917845762
EAN case	4038917845779

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**