



Product description

Order number	70.3040.210
Product description	Pipette tip, working volume: 300 µl, transparent, filling level rings, PCR Performance Tested, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 96 piece(s)/box

Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design
Volume of work	300 µl
Filter	no
Coding colour	orange
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

Max. volume	300 µl
Diameter	7 mm
Length of product	60 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
Certificate of conformity	yes
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order qty.	1920
Type of smallest subpackaging	box
Piece(s) / inner box	480
Piece (s) / outer case	1920
Piece(s) / pallet	61440
Depth of box	530 mm
Width of box	111 mm
Height of box	148 mm
Depth of case	547 mm
Width of case	325 mm
Height of case	245 mm
Case volume	0.0436 cbm
Weight of product	0.0013 kg
Weight of case	4.05 kg
EAN of inner box	4038917433259
EAN case	4038917433242

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

