




## Product Specification

**Filter tip, XL, 1,000 µl, transparent, Biosphere® plus, Low retention, 96 piece(s)/SingleRefill**

**Page 1**



### Product description

Order number	70.3060.375
Product description	Filter tip, working volume: 1,000 µl, transparent, filling level rings, Biosphere® plus, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design, XL, 96 piece(s)/SingleRefill

### Product characteristics

Graduation	yes
Suitable for	SARSTEDT Sarpette® M, Eppendorf, Gilson and Brand and products of identical design
Volume of work	1.000 µl
Filter	yes
Coding colour	blue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

Filter tip, XL, 1,000 µl, transparent, Biosphere® plus, Low retention, 96 piece(s)/SingleRefill

Page 2

### Size

Max. volume	1000 µl
Diameter	9 mm
Length of product	90 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

### Purity & certification

Product category	In vitrodiagnostic device, CE
Certificate of conformity	yes
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

### Packaging

Minimum order qty.	1920
Type of smallest subpackaging	SingleRefill
Piece(s) / inner box	480
Piece (s) / outer case	1920
Piece(s) / pallet	69120
Depth of box	495 mm
Width of box	105 mm
Height of box	148 mm
Depth of case	512 mm
Width of case	325 mm
Height of case	237 mm
Case volume	0.0394 cbm
Weight of product	0.0016 kg
Weight of case	4.54 kg
EAN of inner box	4038917432443
EAN case	4038917432436

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**