



Product Specification
Sarpette® M, single-channel, variable, 100 µl - 1,000 µl, blue, 1 piece(s)/case

Page 1



Product description

Order number	90.3100.000
Product description	Sarpette® M, single-channel, variable, 100 µl - 1,000 µl, blue, 1 piece(s)/case

Product characteristics

Volume of work	100 µl - 1.000 µl
Coding colour	blue
Autoclavable (max. °C)	121 °C

Material & colours

Colour of product	grey
-------------------	------

Purity & certification

Product category	no medical device IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Sarpette® M, single-channel, variable, 100 µl -
1,000 µl, blue, 1 piece(s)/case**

Page 2

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1
Depth of box	315 mm
Width of box	90 mm
Height of box	30 mm
Weight of product	0.088 kg
EAN of inner box	4038917937061
EAN case	4038917937078

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT