



Product description

Order number	90.3108.300
Product description	Sarpette® M, 8-channel, variable, 30 µl - 300 µl, orange, 1 piece(s)/case

Product characteristics

Volume of work	30 µl - 300 µl
Coding colour	orange
Autoclavable (max. °C)	121 °C

Material & colours

Colour of product	grey
-------------------	------

Purity & certification

Product category	no medical device IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1
Depth of box	288 mm
Width of box	150 mm
Height of box	28 mm
Weight of product	0.165 kg
EAN of inner box	4038917937436
EAN case	4038917937443

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

