



### **Product description**

Order number	95.1760.201
Product description	Pipetting aids, storage capacity: 3, for Single and multi-channel pipettes, 1 piece(s)/case

### **Product characteristics**

Suitable for	Single and multi-channel pipettes
Storage capacity	3

### **Material & colours**

Colour of product	white/red
-------------------	-----------

### **Purity & certification**

Product category	no medical device   IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	4
Piece(s) / pallet	144
Depth of box	230 mm
Width of box	150 mm
Height of box	150 mm
Depth of case	360 mm
Width of case	330 mm
Height of case	250 mm
Case volume	0.0297 cbm
Weight of product	0.62 kg
Weight of case	2.92 kg
EAN of inner box	4038917934794
EAN case	4038917934800

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

