




Product Specification
**Serological pipette, with tip, plugged, 1 ml, sterile,
non-pyrogenic/endotoxin-free, non-cytotoxic, 1
piece(s)/blister**



Product description

Order number 86.1251.001
Product description Serological pipette, with tip, plugged, 1 ml, coding colour: yellow, material: PS, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Product characteristics

Label/ Print with printed graduations
Colour of print/label yellow/black
Graduation yes
Design with tip, plugged
Volume of work 0,1 ml - 1,3 ml
Coding colour yellow

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**Serological pipette, with tip, plugged, 1 ml, sterile,
non-pyrogenic/endotoxin-free, non-cytotoxic, 1
piece(s)/blister**

Page 2

Size

Sample volume	1 ml
Diameter	4 mm
Length of product	270 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	blister
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	54000
Depth of case	420 mm
Width of case	250 mm
Height of case	290 mm
Case volume	0.0323 cbm
Weight of product	0.0044 kg
Weight of case	5.01 kg
EAN of inner box	4038917127158
EAN case	4038917021067

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT