




Product Specification

Aspiration pipette, with tip, unplugged, 2 ml,
sterile, non-pyrogenic/endotoxin-free,
non-cytotoxic, 1 piece(s)/blister

Page 1



Product description

Order number	86.1252.011
Product description	Aspiration pipette, with tip, unplugged, 2 ml, material: PS, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Product characteristics

Label/ Print	without print
Graduation	no
Design	with tip, unplugged
Volume of work	0,2 ml - 2 ml

Size

Sample volume	2 ml
Diameter	6 mm
Length of product	270 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**Aspiration pipette, with tip, unplugged, 2 ml,
sterile, non-pyrogenic/endotoxin-free,
non-cytotoxic, 1 piece(s)/blister**

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation Gamma irradiation
Pyrogen-free yes
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest blister
subpackaging
Piece(s) / inner box 100
Piece (s) / outer case 1000
Piece(s) / pallet 54000
Depth of case 420 mm
Width of case 250 mm
Height of case 290 mm
Case volume 0.0323 cbm
Weight of product 0.0044 kg
Weight of case 5.03 kg
EAN of inner box 4038917127202
EAN case 4038917043311

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT