




Product Specification
**Serological pipette, with tip, plugged, 5 ml, sterile,
non-pyrogenic/endotoxin-free, non-cytotoxic, 1
piece(s)/blister**

Page 1



Product description

Order number	86.1253.001
Product description	Serological pipette, with tip, plugged, 5 ml, coding colour: blue, material: PS, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Product characteristics

Label/ Print	with printed graduations
Colour of print/label	blue/black
Graduation	yes
Design	with tip, plugged
Volume of work	0,6 ml - 7,5 ml
Coding colour	blue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**Serological pipette, with tip, plugged, 5 ml, sterile,
non-pyrogenic/endotoxin-free, non-cytotoxic, 1
piece(s)/blister**

Page 2

Size

Sample volume	5 ml
Diameter	8 mm
Length of product	340 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	blister
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	27000
Depth of case	420 mm
Width of case	250 mm
Height of case	290 mm
Case volume	0.0305 cbm
Weight of product	0.008 kg
Weight of case	4.54 kg
EAN of inner box	4038917127240
EAN case	4038917021104

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT