




Product Specification
**Demeter pipette, with tip, plugged, 1.1 ml, sterile,
25 piece(s)/bag**

Page 1



Product description

Order number	86.1686.025
Product description	Demeter pipette, with tip, plugged, 1.1 ml, material: PS, sterile, 25 piece(s)/bag

Product characteristics

Label/ Print	with printed graduations
Colour of print/label	black
Graduation	yes
Design	with tip, plugged
Volume of work	0,5 ml - 1,1 ml

Size

Sample volume	1.1 ml
Diameter	4 mm
Length of product	270 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Demeter pipette, with tip, plugged, 1.1 ml, sterile,
25 piece(s)/bag

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard sterile
Sterilisation Gamma irradiation
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest bag
subpackaging
Piece(s) / inner box 25
Piece (s) / outer case 1000
Piece(s) / pallet 108000
Depth of case 413 mm
Width of case 250 mm
Height of case 150 mm
Case volume 0.0155 cbm
Weight of product 0.0033 kg
Weight of case 3.94 kg
EAN of inner box 4038917133043
EAN case 4038917024662

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT