




Product Specification
**Serological pipette, without tip, plugged, 10 ml,
sterile, 10 piece(s)/bag**

Page 1



Product description

Order number 86.1688.010
Product description Serological pipette, without tip, plugged, 10 ml, material: PS,
sterile, 10 piece(s)/bag

Product characteristics

Label/ Print with printed graduations
Colour of print/label black
Graduation yes
Design without tip, plugged
Volume of work 1 ml - 10 ml

Size

Sample volume 10 ml
Diameter 10 mm
Length of product 335 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Serological pipette, without tip, plugged, 10 ml,
sterile, 10 piece(s)/bag**

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard sterile
Sterilisation Gamma irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of smallest
subpackaging bag
Piece(s) / inner box 10
Piece (s) / outer case 500
Piece(s) / pallet 27000
Depth of case 413 mm
Width of case 250 mm
Height of case 270 mm
Case volume 0.0279 cbm
Weight of product 0.0085 kg
Weight of case 4.97 kg
EAN of inner box 4038917133081
EAN case 4038917024693

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT