




**Product Specification**  
**Transfer pipette, 3.5 ml, (LxW): 155 x 15 mm,**  
**LD-PE, transparent**

**Page 1**



**Product description**

Order number	86.1171.300
Product description	Transfer pipette, nominal volume: 3.5 ml, (LxW): 155 x 15 mm, material: LD-PE, transparent, droplet size: 35 µl - 55 µl, 500 piece(s)/case

**Product characteristics**

Graduation	yes
------------	-----

**Size**

Sample volume	3.5 ml
Width of product	15 mm
Length of product	155 mm
Drop size	35 µl - 55 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Transfer pipette, 3.5 ml, (LxW): 155 x 15 mm,**  
**LD-PE, transparent**

**Page 2**

**Material & colours**

Product material Low Density Polyethylene (LD-PE)  
Colour of product transparent

**Purity & certification**

Product category no medical device | IVD  
Batched yes

**Packaging**

Minimum order qty. 500  
Type of smallest case  
subpackaging  
Piece(s) / inner box 500  
Piece (s) / outer case 5000  
Piece(s) / pallet 56000  
Case volume 0.0994 cbm  
Weight of product 0.0017 kg  
Weight of case 9.41 kg  
EAN of inner box 4038917378383  
EAN case 4038917378376

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**