



Product description

Order number	86.1174.020
Product description	Transfer pipette, nominal volume: 6 ml, (LxW): 152 x 15 mm, material: LD-PE, transparent, droplet size: 40 µl - 55 µl, sterile, 20 piece(s)/bag

Product characteristics

Graduation	no
------------	----

Size

Sample volume	6 ml
Width of product	15 mm
Length of product	152 mm
Drop size	40 µl - 55 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Low Density Polyethylene (LD-PE)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 3200
Type of smallest bag
subpackaging
Piece(s) / inner box 400
Piece (s) / outer case 3200
Piece(s) / pallet 48000
Case volume 0.1395 cbm
Weight of product 0.0023 kg
Weight of case 8.95 kg
EAN of inner box 4038917318402
EAN case 4038917318396

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

