



Product description

Order number	86.1174.300
Product description	Transfer pipette, nominal volume: 6 ml, (LxW): 152 x 15 mm, material: LD-PE, transparent, droplet size: 40 μ l - 55 μ l, 400 piece(s)/case

Product characteristics

Graduation	no
------------	----

Size

Sample volume	6 ml
Width of product	15 mm
Length of product	152 mm
Drop size	40 μ l - 55 μ l

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Transfer pipette, 6 ml, (LxW): 152 x 15 mm, LD-PE,
transparent**

Page 2

Material & colours

Product material Low Density Polyethylene (LD-PE)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 400
Type of smallest case
subpackaging
Piece(s) / inner box 400
Piece (s) / outer case 4000
Piece(s) / pallet 44800
Case volume 0.0994 cbm
Weight of product 0.002 kg
Weight of case 9.07 kg
EAN of inner box 4038917378505
EAN case 4038917378499

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT