



Product description

Order number	86.1175.020
Product description	Transfer pipette, nominal volume: 6 ml, (LxW): 146 x 15 mm, material: LD-PE, transparent, droplet size: 20 μ l - 30 μ l, sterile, 20 piece(s)/bag

Product characteristics

Graduation no

Size

Sample volume	6 ml
Width of product	15 mm
Length of product	146 mm
Drop size	20 μ l - 30 μ l

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Low Density Polyethylene (LD-PE)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 3200
Type of smallest subpackaging bag
Piece(s) / inner box 400
Piece (s) / outer case 3200
Piece(s) / pallet 3200
Case volume 0.1309 cbm
Weight of product 0.0022 kg
Weight of case 8.64 kg
EAN of inner box 4038917318433
EAN case 4038917318426

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

