



Product description

Order number 86.1568.050
Product description Inoculation needle, material: PS, orange, sterile, 50 piece(s)/blister

Size

Width of product 4 mm
Length of product 197 mm

Material & colours

Product material Polystyrene (PS)
Colour of product orange

Purity & certification

Product category no medical device | IVD
Purity standard sterile
Sterilisation Gamma irradiation
Batched yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	2000
Type of smallest subpackaging	blister
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	264000
Depth of box	216 mm
Width of box	161 mm
Height of box	101 mm
Depth of case	338 mm
Width of case	228 mm
Height of case	120 mm
Case volume	0.0092 cbm
Weight of product	0.0012 kg
Weight of case	2.79 kg
EAN of inner box	4038917130974
EAN case	4038917023542

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

