



Product description

Order number	82.1195
Product description	Petri dish, (ØxH): 92 x 16 mm, material: PS, transparent, 2 compartments, with ventilation cams, 20 piece(s)/tube bag

Product characteristics

Minimum temperature	0 °C
Maximum temperature	80 °C
Ventilation cams	yes

Size

Max. volume	60 ml
Diameter	92 mm
Height of product	16 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Petri dish, 92 x 16 mm, transparent, 2 compartments, with ventilation cams

Page 2

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Latex-free yes
Safety data sheet no
Batched yes

Packaging

Minimum order qty. 480
Type of smallest subpackaging tube bag
Piece(s) / inner box 480
Piece (s) / outer case 480
Piece(s) / pallet 9600
Depth of case 584 mm
Width of case 387 mm
Height of case 360 mm
Case volume 0.0814 cbm
Weight of product 0.0163 kg
Weight of case 8.93 kg
EAN of inner box 4038917085700
EAN case 4038917004336

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT