



### **Product description**

Order number	82.1473.040
Product description	Petri dish, (ØxH): 92 x 16 mm, material: PS, yellow, with ventilation cams, 20 piece(s)/tube bag

### **Product characteristics**

Minimum temperature	0 °C
Maximum temperature	80 °C
Ventilation cams	yes

### **Size**

Max. volume	80 ml
Diameter	92 mm
Height of product	16 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material Polystyrene (PS)  
Colour of product yellow

**Purity & certification**

Product category no medical device | IVD  
Latex-free yes  
Safety data sheet no  
Batched yes

**Packaging**

Minimum order qty. 480  
Type of smallest tube bag  
subpackaging  
Piece(s) / inner box 480  
Piece (s) / outer case 480  
Piece(s) / pallet 9600  
Depth of case 584 mm  
Width of case 387 mm  
Height of case 360 mm  
Case volume 0.0814 cbm  
Weight of product 0.0169 kg  
Weight of case 9.25 kg  
EAN of inner box 4038917129374  
EAN case 4038917022156

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

