



Product description

| | |
|---------------------|--|
| Order number | 82.9923.426 |
| Product description | Petri dish, (ØxH): 54.65 x 14.7 mm, material: PS, transparent, without ventilation cams, 10 piece(s)/bag |

Product characteristics

| | |
|------------------|----|
| Ventilation cams | no |
|------------------|----|

Size

| | |
|-------------------|----------|
| Diameter | 54.65 mm |
| Height of product | 14.7 mm |

Material & colours

| | |
|-------------------|------------------|
| Product material | Polystyrene (PS) |
| Colour of product | transparent |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

| | |
|-------------------|-----|
| Latex-free | yes |
| Pyrogen-free | no |
| Safety data sheet | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 10 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 40000 |
| Depth of case | 604 mm |
| Width of case | 264 mm |
| Height of case | 373 mm |
| Case volume | 0.06 cbm |
| Weight of case | 8.2 kg |
| EAN of inner box | 4038917154208 |
| EAN case | 4038917039772 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

