



### **Product description**

|                     |   |
|---------------------|---|
| Order number        | 82.9923.429   |
| Product description | Petri dish, (ØxH): 68.45 x 15 mm, material: PS, transparent, without ventilation cams, sterile, 10 piece(s)/bag |

### **Product characteristics**

|                  |    |
|------------------|----|
| Ventilation cams | no |
|------------------|----|

### **Size**

|                   |          |
|-------------------|----------|
| Diameter          | 68.45 mm |
| Height of product | 15 mm    |

### **Material & colours**

|                   |                  |
|-------------------|------------------|
| Product material  | Polystyrene (PS) |
| Colour of product | transparent      |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

|                   |         |
|-------------------|---------|
| Latex-free        | yes     |
| Purity standard   | sterile |
| Pyrogen-free      | no      |
| Safety data sheet | no      |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 600           |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 10            |
| Piece (s) / outer case        | 600           |
| Piece(s) / pallet             | 28800         |
| Depth of case                 | 608 mm        |
| Width of case                 | 264 mm        |
| Height of case                | 324 mm        |
| Case volume                   | 0.0519 cbm    |
| Weight of case                | 6.7 kg        |
| EAN of inner box              | 4038917154239 |
| EAN case                      | 4038917039802 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

